

# **Cambridge International Examinations**

Cambridge International General Certificate of Secondary Education

BIOLOGY 0610/52

Paper 5 Practical Test

October/November 2015

CONFIDENTIAL INSTRUCTIONS

Great care should be taken to ensure that any confidential information given does not reach the candidates either directly or indirectly.

If you have any problems or queries regarding these Confidential Instructions, please contact CIE

by e-mail info@cie.org.uk by phone +44 1223 553554 by fax +44 1223 553558

stating the Centre number, the nature of the query and the syllabus number quoted above.



#### READ THESE INSTRUCTIONS FIRST

These instructions give details of the apparatus required by each candidate for each experiment in this paper. A summary of the questions that will be presented to the candidates is included, where appropriate, to allow the biology teacher to test the apparatus appropriately. **Testing must be done out of sight of all candidates.** 

No access to the question paper is permitted in advance of the examination.

Centres are reminded that candidates are expected to follow the instructions on the question paper and record all their results. They will not be penalised if these results are not what they expect.

The Supervisor should make sure that the Supervisor's Report is fully completed and a copy is enclosed with each packet of scripts.

It is assumed that the ordinary apparatus of a science laboratory will be available, including a supply of purified water (distilled or deionised).

If arrangements are made for different sessions for different groups of candidates, care must be taken to ensure that the different groups of candidates are effectively isolated so that **no information passes** between them.

All specimens should carry only the code letters and numbers as indicated and their identity should not be revealed to the candidates.

Supervisors should ensure that all specimens have the correct identity attached to the specimen and that these are **not** removed during the examination.

If a candidate breaks any of the apparatus, or loses any of the material supplied, the matter should be rectified and a note made in the Supervisor's Report.

Supervisors are advised to remind candidates that **all** substances in the examination should be treated with caution. Pipette fillers and safety goggles should be used where necessary.

In accordance with the COSHH (Control of Substances Hazardous to Health) Regulations, operative in the UK, a hazard appraisal of the examination has been carried out.

The following codes are used where relevant:

 $\mathbf{C}$  = corrosive substance  $\mathbf{F}$  = highly flammable substance

**H** = harmful or irritating substance **O** = oxidising substance

N =harmful to the environment T =toxic substance

Centres are reminded that they are **not** permitted to open the question paper envelopes before the examination.

If there are any difficulties with any aspect of setting up this practical examination that the Centre is not able to resolve, it is essential for Centres to contact the Product Manager as soon as possible by **e-mail** to info@cie.org.uk, by **fax** to +44 1223 553558 or by **phone** to +44 1223 553554.

#### Question 1

Each candidate should be provided with:

- (i) three large test-tubes (approximately  $25 \,\mathrm{mm} \times 150 \,\mathrm{mm}$ )
- (ii) three standard test-tubes (approximately 15 mm × 150 mm)
- (iii) means of measuring 10 cm<sup>3</sup>, e.g. 10 cm<sup>3</sup> syringe
- (iv) two  $\times$  means of measuring 1 cm<sup>3</sup>, e.g. two 1 cm<sup>3</sup> syringes
- (v) one large empty beaker (capable of supporting three large and three standard test-tubes)
- (vi) three beakers, each deep enough to support a large test-tube and a standard test-tube:
  - a. one containing an ice and water mixture, labelled cold
  - b. one containing water at room temperature, labelled warm
  - c. one empty beaker labelled hot.

Candidates have been asked to raise their hand when they need hot water.

- (vii) supply of hot water at 55°C
- (viii) one thermometer, range –10 °C to 100 °C
  - (ix) 10 cm<sup>3</sup> yeast suspension in a closed container, labelled yeast suspension
  - (x) 10 cm<sup>3</sup> of 0.005 % methylene blue solution, labelled methylene blue solution
  - (xi) 40 cm<sup>3</sup> of 10% glucose solution, labelled glucose solution
- (xii) glass rod
- (xiii) access to paper towels (in case of spillages)
- (xiv) view of a clock or timer that allows the candidates to time seconds accurately.

#### **Yeast stock suspension preparation** (sufficient for 10 candidates):

- Weigh 20 g of dry active baking yeast. This should be new stock.
- Add 100 cm<sup>3</sup> warm water\*. This should be at approximately 35°C.
- Stir thoroughly.

Sufficient solution for the candidates should be made about 15-30 minutes before the start of the examination and left to start to respire. This can be kept at room temperature once mixed.

Please note that there is great variation between yeast sold for baking purposes. Brewer's yeast is NOT suitable. It is essential that the yeast is active and respiring, so trialling must be carried out in sufficient time to allow sourcing of a suitable yeast (see page 4 for details of trialling).

<sup>\*</sup>If tap water is used, it should be boiled before use.

#### Solution preparation:

# 0.1% methylene blue stock solution

- Weigh 0.1 g methylene blue-trihydrate (*M*, 373.9).
- Add 50 cm<sup>3</sup> of distilled water and stir until dissolved.
- Make up to 100 cm<sup>3</sup> with distilled water.

This should be prepared a week before the examination and stored in a cool place.

## **0.005% methylene blue** (sufficient for 20 candidates)

- Take 10 cm<sup>3</sup> of 0.1% methylene blue stock solution and make up to 200 cm<sup>3</sup> with distilled water.
- Stir to mix.

This should be made the day before the examination and stored in a cool place.

**Before the examination**, the Supervisor must try out this part of the investigation.

The Supervisor should carry out part of the investigation by adding  $10\,\mathrm{cm^3}$  of 10% glucose solution and  $1\,\mathrm{cm^3}$  of 0.005% methylene blue solution to a large test-tube.  $1\,\mathrm{cm^3}$  of the yeast suspension should be added to a standard test-tube. Both tubes should be kept at room temperature ( $20^\circ\mathrm{C}-25^\circ\mathrm{C}$ ) for 5 minutes. The yeast suspension should then be tipped into the large test-tube containing the mixture of glucose and methylene blue solutions and shaken until the blue colour is evenly spread. The time for the methylene blue colour to fade completely should be between 10-12 minutes at room temperature. If a longer time is recorded with yeast that is actively respiring, then the yeast concentration can be increased by up to a further 5%.

It may be necessary to try more than one type of dried yeast as the activity can vary greatly.

The Supervisor (**not** the Invigilator) should carry out the practical aspects of **Question 1** and record their results in the space in the Supervisor's Report. This must be done during the examination, using the same apparatus and reagents as the candidates but **out of sight** of the candidates.

### **Question 2**

Each candidate should be provided with:

(i) one ruler (with mm scale)

# **BLANK PAGE**

6

### **BLANK PAGE**

Permission to reproduce items where third-party owned material protected by copyright is included has been sought and cleared where possible. Every reasonable effort has been made by the publisher (UCLES) to trace copyright holders, but if any items requiring clearance have unwittingly been included, the publisher will be pleased to make amends at the earliest possible opportunity.

To avoid the issue of disclosure of answer-related information to candidates, all copyright acknowledgements are reproduced online in the Cambridge International Examinations Copyright Acknowledgements Booklet. This is produced for each series of examinations and is freely available to download at www.cie.org.uk after the live examination series.

Cambridge International Examinations is part of the Cambridge Assessment Group. Cambridge Assessment is the brand name of University of Cambridge Local Examinations Syndicate (UCLES), which is itself a department of the University of Cambridge.

# This form should be completed and returned to CIE with the scripts.

# **REPORT ON PRACTICAL BIOLOGY**

The Supervisor should provide the following information.

1	Was	any difficulty experienced in providing the necessary materials? If so, give brief details.
2		e details of any difficulties experienced by particular candidates, giving names and candidate obers. Reference should be made to:
	(a)	difficulties with specimens or materials;
	(b)	accidents to apparatus or materials;
	(c)	assistance provided in case of colour-blindness;
	(d)	any other information that is likely to assist the Examiner, especially if this cannot be discovered from the scripts.
		er cases of individual hardship, e.g. illness or disability, should be reported direct to CIE on the nal 'Special Consideration Form' as detailed in the Handbook for Centres.
3	prac Res	ing the examination, the Supervisor should, <b>out of sight of the candidates</b> , carry out the stical aspects of <b>Question 1</b> , using the same apparatus and reagents as the candidates. ults should be recorded in the space on page 8 ( <b>not</b> on a spare Question Paper). Attach extra ets if necessary.
	The	Invigilator should <b>not</b> carry out <b>Question 1</b> .

	plan of work benches, giving details of the candidate numbers of the places occupied by the andidates for each session, must be enclosed with the scripts.
Declar	ration (to be signed by the Principal or the Examinations Officer)
The present	reparation of this practical has been carried out so as to maintain fully the security of the nation.
	Signed
	Name (in block capitals)
	Centre number
Centre	e name
	ots are required to be dispatched in more than one packet, it is essential that a copy of the visor's Report and the appropriate seating plan(s) are inside <b>each packet</b> .

© UCLES 2015 0610/52/CI/O/N/15